

RIH – ABDOMEN PELVIS ANGIO FOR LOWER GI BLEED GE LIGHTSPEED VCT PROTOCOL

Indications: For rapid detection of lower gastrointestinal bleeding.

| | | | | |
|--|---|------------------------------------|-----------------------|----------------------|
| Position/Landmark | Head first or feet first-Supine Xyphoid | | | |
| Topogram Direction | Craniocaudal | | | |
| Respiratory Phase | Inspiration | | | |
| Scan Type | Helical | | | |
| KV / mA / Rotation time (sec) Pitch / Speed (mm/rotation) Noise Index / ASiR / Dose Reduction | 120kv / smart mA (120-450) / 0.5 sec .984:1 , 39.37mm 11.5 nc and 16 contrast / 70 / 30% | | | |
| Detector width x Rows = Beam Collimation | 0.625mm x 64 = 40mm | | | |
| Average Tube Output | Each Helical: ctdi – 11.3mGy dlp – 616 mGy.cm | | | |
| First Helical Set Slice Thickness/ Spacing Algorithm Recon Destination | recon | body part | thickness/ spacing | recon destination . |
| | 1 | nc abd/pelvis | 5mm x 5mm | standard pacs |
| Second Helical Set Slice Thickness/ Spacing Algorithm Recon Destination | recon | body part | thickness/ spacing | recon destination . |
| | 1 | abd pelvis ct angio | 2.5mm x 2.5mm | standard pacs |
| | 2 | thin ct angio | .6mm x .6mm | standard for dmpr |
| Third Helical Set Slice Thickness/ Spacing Algorithm Recon Destination | recon | body part | thickness/ spacing | recon destination . |
| | 1 | delayed abd pelvis | 2.5mm x 2.5mm | standard pacs |
| | 2 | thin delayed ct | .6mm x .6mm | standard for dmpr |
| Scan Start / End Locations DFOV | 1 cm superior to diaphragm lesser trochanters 38cm decrease appropriately | | | |
| IV Contrast Volume / Type / Rate | 100mL Iohexol (Omnipaque 350) 4mL/sec | | | |
| Scan Delay | Non-Contrast ----- | CTA smart prep at celiac artery | Delay 80 seconds | |
| 2D/3D Technique Used | CTA: DMPR of 2.5 mm x 2.5 mm coronal abdomen/pelvis series (auto-batch on), mip mode . 10mm x 1mm inverted coronal abdomen/pelvis mip series (auto-batch on) Delay: DMPR of 2.5 mm x 2.5 mm coronal abdomen/pelvis series (auto-batch on), mip mode , auto-transferred to PACS | | | |
| Comments: | A non-contrast study is done first. Then the cta is done using a smart prep at the level of the celiac artery. Note: There is a second helical scan done 60 seconds after the cta to look for subtle blood pooling. The arterial series has an inverted coronal mip series. The gray scales should be inverted in pacs for this series only. | | | |
| Images required in PACS | Scouts, 5mm x 5mm axial nc abdomen/pelvis, 2.5mm x 2.5mm axial cta abdomen/pelvis, 2.5mm x 2.5mm coronal mip arterial abdomen/pelvis, 10mm x 1mm inverted coronal abdomen/pelvis mip, 2.5mm x 2.5mm axial delayed abdomen/pelvis, 2.5mm x 2.5mm coronal mip delayed abdomen/pelvis, Dose Report | | | |