

# RIH – LUNG SCREENING GE LIGHTSPEED VCT PROTOCOL

**Indications – Survey of the lungs for nodules.**

<b>Position/Landmark</b>	Head first or feet first-Supine Sternal Notch				
<b>Respiratory Phase</b>	Inspiration				
<b>Topogram Direction</b>	Craniocaudal				
<b>Scan Type</b>	Helical				
<b>KV / mA / Rotation time (sec) Pitch / Speed (mm/rotation) Noise Index / ASiR / Dose Reduction</b>	120kv / smart mA (50-150) / 0.5 sec 1.375:1 , 27.5mm 33.0 / 70 / 0%				
<b>Detector width x Rows = Beam Collimation</b>	0.625mm x 32 = 20mm				
<b>Average Tube Output</b>	ctdi – 4 mGy dlp – 110 mGy.cm				
<b>Helical Set</b>					recon
Slice Thickness/ Spacing	recon	body part	thickness/ spacing	algorithm	recon destination .
Algorithm	1	<b>chest</b>	2.5mm x 2.5mm	standard	pacs
Recon Destination	2	thin chest	.6mm x .6mm	standard	for dmpr
	3	<b>lung</b>	2.5mm x 2.5mm	lung	pacs
<b>Scan Start / End Locations</b>	1cm superior to lung apices 1cm inferior to costophrenic angles				
<b>DFOV</b>	38cm decrease appropriately				
<b>IV Contrast Volume / Type / Rate</b>					
<b>Scan Delay</b>					
<b>2D/3D Technique Used</b>	DMPR of 2.5mm x 2.5mm <b>coronal chest</b> series (auto-batch on), average mode, auto-transferred to PACS.				
<b>Comments:</b> Recon 2 is a single thin helical group of the chest for direct mpr.					
<b>Images required in PACS</b>	Scouts, 2.5mm x 2.5mm axial chest, 2.5mm x 2.5mm coronal chest, 2.5mm x 2.5mm axial lungs, Dose Report				