

RIH – PELVIS FOR FRACTURE GE LIGHTSPEED VCT PROTOCOL

Indication: trauma, fracture, dislocation, abscess

Position/Landmark	Head first or feet first-Supine Iliac Crest			
Topogram Direction	Craniocaudal			
Respiratory Phase	Suspension			
Scan Type	Helical			
KV / mA / Rotation time (sec)	120kv / smart mA (120-450) / 0.5 sec			
Pitch / Speed (mm/rotation)	.984:1 , 39.37mm			
Noise Index / ASiR / Dose Reduction	11.5 / 70 / 30%			
Detector width x Rows = Beam Collimation	0.625mm x 64 = 40mm			
Average Tube Output	ctdi – 10.7 mGy dlp – 313 mGy.cm			
Helical Set	recon	body part	thickness/ spacing	recon destination .
Slice Thickness/ Spacing			algorithm	
Algorithm	1	pelvis bone	2.5mm x 2.5mm	bone pacs
Recon Destination	2	thin pelvis/hip	.6mm x .6mm	bone for dmpr
	3	pelvis soft tissue	5mm x 5mm	standard pacs
Scan Start / End Locations	1 cm superior to iliac crest lesser trochanters			
DFOV	38cm decrease appropriately			
IV Contrast Volume / Type / Rate				
Scan Delay				
2D/3D Technique Used	DMPR: sagittal and coronal reformats , 3.0mm x 3.0mm, average mode			
Comments: Recon 2 is a thin helical volume of the pelvis that is archived and used in direct multi-planar reformats.				
When a ct cystogram is ordered, instill 50mL of Iohexol (Omnipaque 240) into 500mL of normal saline and retrograde drip into the bladder via the patient's foley.				
Images required in PACS	Scouts, 2.5mm x 2.5mm axial pelvis for bone, 5mm x 5mm pelvis for soft tissue, 3mm x 3mm sagittal reformats, 3mm x 3mm coronal reformats, Dose Report			